

cosine[®]45

艾柯医疗
ACCUMEDICAL

远端通路导管



cosine[®]45

远端通路导管



更广
兼容性



更高
到位性



跨迂曲
能力



可抽吸



多术式
场景

规格

对应特性



外径 (通体均一)

0.060" | 1.52mm | 4.6F

可在 $\geq 6F$ 通路中穿行:

6F短鞘、长鞘, 6F普通导管、中间导管




内径 (通体均一)

0.046" | 1.17mm | 3.5F

内部可通行所有单导管系统:

0.017" / 0.021" / 0.027" / 0.029" 微导管



长度

70~140cm

可到达更远的血管节段:


M2/3, A2/3, P2/3

型号规格	有效长度(cm)	型号规格	有效长度(cm)	内径	外径
SC401SHR	70	SC413SHR	118	0.046" 1.17mm 3.5F	0.060" 1.52mm 4.6F
SC402SHR	75	SC414SHR	120		
SC403SHR	80	SC415SHR	122		
SC404SHR	85	SC416SHR	124		
SC405SHR	90	SC417SHR	126		
SC406SHR	95	SC418SHR	128		
SC407SHR	100	SC419SHR	130		
SC408SHR	105	SC420SHR	132		
SC409SHR	110	SC421SHR	134		
SC410SHR	112	SC422SHR	136		
SC411SHR	114	SC423SHR	138		
SC412SHR	116	SC424SHR	140		

应用场景


1

极远端/迂曲路径的单纯
弹簧圈栓塞术




2

迂曲路径/宽颈巨大动脉瘤的
FD栓塞术




3

WEB栓塞术




4

MeVO(中型动脉栓塞)
机械取栓术



5

颅内狭窄
球囊/支架成形术



艾柯医疗器械（北京）股份有限公司
北京市朝阳区东三环工人体育场北路甲2号盈科中心A栋1210



www.accumedical.com
+86 010 6506 5018
