

# sioe21



微导管

艾柯医疗  
ACCUMEDICAL



# sioe21微导管

详细规格



A: 内径 = 0.021", 0.53mm  
 B: 长度 = 100~154cm  
 C: 外径 = 0.034", 0.86mm  
 D: 亲水涂层长度 = 95cm

## 参数配置

国械注20233030768

型号规格	内径	外径	有效长度
MC2101HR	0.021" 0.53mm	0.034" 0.86mm	100cm
MC2102HR			105cm
MC2103HR			110cm
MC2104HR			115cm
MC2105HR			120cm
MC2106HR			125cm
MC2107HR			130cm
MC2108HR			135cm
MC2109HR			140cm
MC2110HR			145cm
MC2111HR			150cm
MC2112HR			151cm
MC2113HR			152cm
MC2114HR			153cm
MC2115HR			154cm