


cosineTM

艾柯医疗
ACCUMEDICAL

远端通路导管




顺滑到位 稳行致远

国械注准20223030188

cosine® 远端通路导管

详细规格



A: 内径 = 0.055", 1.4mm / 0.071", 1.8mm
B: 长度 = 70~138cm
C: 外径 = 0.068", 1.73mm / 0.084", 2.13mm
D: 亲水涂层长度 = 40cm

参数配置

国械注准20223030188

型号规格		外径		内径		有效长度
支撑型 5F	SC502SHR	0.068"	1.73mm	0.055"	1.4mm	75cm
	SC504SHR	0.068"	1.73mm	0.055"	1.4mm	85cm
	SC506SHR	0.068"	1.73mm	0.055"	1.4mm	95cm
	SC508SHR	0.068"	1.73mm	0.055"	1.4mm	105cm
	SC510SHR	0.068"	1.73mm	0.055"	1.4mm	115cm
	SC516SHR	0.068"	1.73mm	0.055"	1.4mm	125cm
	SC521SHR	0.068"	1.73mm	0.055"	1.4mm	130cm
	SC526SHR	0.068"	1.73mm	0.055"	1.4mm	135cm
	SC529SHR	0.068"	1.73mm	0.055"	1.4mm	138cm
支撑型 6F	SC601SHR	0.084"	2.13mm	0.071"	1.8mm	70cm
	SC603SHR	0.084"	2.13mm	0.071"	1.8mm	80cm
	SC605SHR	0.084"	2.13mm	0.071"	1.8mm	90cm
	SC608SHR	0.084"	2.13mm	0.071"	1.8mm	105cm
	SC610SHR	0.084"	2.13mm	0.071"	1.8mm	115cm
	SC611SHR	0.084"	2.13mm	0.071"	1.8mm	120cm
	SC616SHR	0.084"	2.13mm	0.071"	1.8mm	125cm
	SC621SHR	0.084"	2.13mm	0.071"	1.8mm	130cm
	SC626SHR	0.084"	2.13mm	0.071"	1.8mm	135cm
柔软型 5F	SC501SHX	0.068"	1.73mm	0.055"	1.4mm	70cm
	SC503SHX	0.068"	1.73mm	0.055"	1.4mm	80cm
	SC505SHX	0.068"	1.73mm	0.055"	1.4mm	90cm
	SC506SHX	0.068"	1.73mm	0.055"	1.4mm	95cm
	SC508SHX	0.068"	1.73mm	0.055"	1.4mm	105cm
	SC510SHX	0.068"	1.73mm	0.055"	1.4mm	115cm
	SC511SHX	0.068"	1.73mm	0.055"	1.4mm	120cm
	SC516SHX	0.068"	1.73mm	0.055"	1.4mm	125cm
	SC521SHX	0.068"	1.73mm	0.055"	1.4mm	130cm
柔软型 6F	SC601SHX	0.084"	2.13mm	0.071"	1.8mm	70cm
	SC603SHX	0.084"	2.13mm	0.071"	1.8mm	80cm
	SC605SHX	0.084"	2.13mm	0.071"	1.8mm	90cm
	SC607SHX	0.084"	2.13mm	0.071"	1.8mm	100cm
	SC610SHX	0.084"	2.13mm	0.071"	1.8mm	115cm
	SC616SHX	0.084"	2.13mm	0.071"	1.8mm	125cm
	SC626SHX	0.084"	2.13mm	0.071"	1.8mm	135cm
	SC629SHX	0.084"	2.13mm	0.071"	1.8mm	138cm