



attractor[®]

Aspiration Catheter



attractor[®] Aspiration Catheter

Specifications



A: Inner Diameter (ID) = 0.071in, 1.8mm
 B: Length = 90–138cm
 C: Outer Diameter (OD) = 0.084in, 2.13mm, 6F
 D: Hydrophilic Coating Length = 40cm

Medical Device Registration Certificate No. of P.R. China: 20233031587

SKU	OD		ID		Effective Length
	Proximal end	Distal end	Proximal end	Distal end	
ACP01STR					90cm
ACP02STR					95cm
ACP03STR					100cm
ACP04STR					105cm
ACP05STR					110cm
ACP06STR					115cm
ACP07STR					120cm
ACP08STR					121cm
ACP09STR					122cm
ACP10STR					123cm
ACP11STR					124cm
ACP12STR					125cm
ACP13STR	2.13mm	2.13mm	1.80mm	1.80mm	126cm
ACP14STR	0.084in	0.084in	0.071in	0.071in	127cm
ACP15STR					128cm
ACP16STR					129cm
ACP17STR					130cm
ACP18STR					131cm
ACP19STR					132cm
ACP20STR					133cm
ACP21STR					134cm
ACP22STR					135cm
ACP23STR					136cm
ACP24STR					137cm
ACP25STR					138cm