

cosine[®]45

艾柯医疗
ACCUMEDICAL

Distal Access Catheter



cosine[®]45

Distal Access Catheter



Compatibility



Distal Access



Aspiration



Trackable



Various Applications

Specifications

Characteristics



OD (uniform)
0.060" | 1.52mm | 4.6F

≥6F Guiding:
6F short sheath and long sheath; 6F guiding catheter



ID (uniform)
0.046" | 1.17mm | 3.5F

Compatible with all microcatheters:
0.017" / 0.021" / 0.027" / 0.029"



Length
70~140cm


More distal access:
M2/3, A2/3, P2/3

SKU	Effective Length (cm)	SKU	Effective Length (cm)	Inner Diameter	Outer Diameter
SC401SHR	70	SC413SHR	118	0.046" 1.17mm 3.5F	0.060" 1.52mm 4.6F
SC402SHR	75	SC414SHR	120		
SC403SHR	80	SC415SHR	122		
SC404SHR	85	SC416SHR	124		
SC405SHR	90	SC417SHR	126		
SC406SHR	95	SC418SHR	128		
SC407SHR	100	SC419SHR	130		
SC408SHR	105	SC420SHR	132		
SC409SHR	110	SC421SHR	134		
SC410SHR	112	SC422SHR	136		
SC411SHR	114	SC423SHR	138		
SC412SHR	116	SC424SHR	140		

Various Applications


1

Coiling for Aneurysms in Distal or Tortuous Vessels




2

Flow Diverter for Aneurysms in Distal or Tortuous Vessels




3

WEB for Aneurysms




4

Mechanical Thrombectomy in MeVO (Medium Vessel Occlusions)



5

Balloon Angioplasty / Stenting for Intracranial Stenosis



AccuMedical Beijing Inc.
1210, Building A Pacific Century Place, 2A Gong Ti Bei Lu, Chaoyang
District, Beijing, CHINA



www.accumedical.com
+86 010 6506 5018
