



sioe 21

Microcatheter



sioe21 Microcatheter

Specifications



A: Inner Diameter (ID) = 0.021", 0.53mm
 B: Length = 100~154cm
 C: Outer Diameter (OD) = 0.034", 0.86mm
 D: Hydrophilic Coating Length = 95cm

Medical Device Registration Certificate No. of P.R. China: 20233030768

SKU	ID	OD	Effective Length
MC2101HR			100cm
MC2102HR			105cm
MC2103HR			110cm
MC2104HR			115cm
MC2105HR			120cm
MC2106HR			125cm
MC2107HR			130cm
MC2108HR	0.021"	0.034"	135cm
MC2109HR	0.53mm	0.86mm	140cm
MC2110HR			145cm
MC2111HR			150cm
MC2112HR			151cm
MC2113HR			152cm
MC2114HR			153cm
MC2115HR			154cm