

sioe™


微导管

艾柯医疗
ACCUMEDICAL



推送无阻 弦无虚发

国械注准20223030811



A: 内径 = 0.027", 0.69mm
B: 长度 = 90~154cm
C: 外径 = 0.041", 1.02mm
D: 亲水涂层长度 = 85cm

参数配置

国械注准20223030811

型号规格	外径		内径		有效长度
MC2705HR	3.1F	1.02mm	0.027"	0.69mm	90cm
MC2706HR	3.1F	1.02mm	0.027"	0.69mm	95cm
MC2707HR	3.1F	1.02mm	0.027"	0.69mm	100cm
MC2708HR	3.1F	1.02mm	0.027"	0.69mm	105cm
MC2709HR	3.1F	1.02mm	0.027"	0.69mm	110cm
MC2710HR	3.1F	1.02mm	0.027"	0.69mm	115cm
MC2711HR	3.1F	1.02mm	0.027"	0.69mm	120cm
MC2712HR	3.1F	1.02mm	0.027"	0.69mm	125cm
MC2713HR	3.1F	1.02mm	0.027"	0.69mm	130cm
MC2714HR	3.1F	1.02mm	0.027"	0.69mm	135cm
MC2715HR	3.1F	1.02mm	0.027"	0.69mm	140cm
MC2717HR	3.1F	1.02mm	0.027"	0.69mm	150cm
MC2721HR	3.1F	1.02mm	0.027"	0.69mm	154cm